

# ELEVATE<sup>TM</sup>MS

Take TMS to the next level

For 20 years, Rogue Research has worked to develop tools that help you advance the boundaries of non-invasive brain stimulation.

We expand this tradition of offering the best tool possible into the clinical domain with the release of our latest product, **ELEVATE TMS**.


Based on a completely different architecture than traditional TMS, called cTMS, **ELEVATE TMS** offers more **control** than any TMS device available today.



**Rogue Research Inc.**

## Features of **ELEVATE**TMS

- First new TMS design in years
- Variable pulse shapes to reliably excite or inhibit neuronal circuits
- Monophasic, biphasic, unidirectional pulses
- Integrated output recording
- External control via sequence file and Ethernet interface



## Pulse Waveform Capabilities

- Directionality Control: Controls amplitude of the -ve phase relative to the +ve phase
- Variable pulse widths up to 385  $\mu$ Sec (with optional high inductance coil)
- Monophasic, biphasic, polyphasic, staircase, asymmetric
- Repetition rates up to 1 kHz
- Unidirectional theta burst
- Charging power: 2x1500 Joules per second

## One of the Key Components in NIBS



### Targeting

Use the most appropriate method to determine the target location for stimulation according to your research goals.



### Dosing

**Elevate TMS** gives you control over how the TMS pulse modulates cortical excitability. Select the best pulse shape to reliably achieve the desired effect.



### Execution

Use any method to place the coil over the location, including neuronavigation (not included).

**ELEVATE**TMS




## Modern, intuitive user interface

- Large, capacitive touch screen
- Simple and advanced controls
- Predicted waveform displayed initially and actual measured output overlaid for confirmation of delivered pulse



## Thoughtful Coil Designs

- Removable handles to suit every preference
- Different inductances to extend cTMS pulse range
- B-Field 3D mapped for use in E-field modelling research
- Integrated coil tracker mount
- More models to come (cooled coil, 50mm fig-8 etc.)



The Journal of Physiology 596(20) June 2018 DOI: 10.1113/JP275798

**Effect of Coil Orientation on Motor Evoked Potentials in Humans with Tetraplegia**

Hang Jin JoHang Jin JoVincenzo Di LazzaroMonica A. Perez

Brain Stimul. 2016 Jan-Feb;9(1):39-47. doi: 10.1016/j.brs.2015.08.013. Epub 2015 Sep 1.

**Enhancement of Neuromodulation with Novel Pulse Shapes Generated by Controllable Pulse Parameter Transcranial Magnetic Stimulation.**

Goetz SM1, Luber B2, Lisanby SH2, Murphy DL1, Kozyrkov IC1, Grill WM3, Peterchev AV4.

Clin Neurophysiol, 127(1), 675-83 (2015)

**Effect of coil orientation on strength-duration time constant and I-wave activation with controllable pulse parameter transcranial magnetic stimulation.**

D'Ostilio K, Goetz SM, Hannah R, Ciocca M, Chieffo R, Chen JC, Peterchev AV, Rothwell JC

Clin Neurophysiol. 2016 Jan;127(1):675-83. doi: 10.1016/j.clinph.2015.05.017. Epub 2015 May 30.

**Effect of coil orientation on strength-duration time constant and I-wave activation with controllable pulse parameter transcranial magnetic stimulation.**

D'Ostilio K1, Goetz SM2, Hannah R3, Ciocca M4, Chieffo R5, Chen JC6, Peterchev AV7, Rothwell JC3.

Clin Neurophysiol 125, S1-S339Peterchev, A.V., D'Ostilio, K., Rothwell, J.C., Murphy, D.L. (2014).

**Intermittent theta burst stimulation inhibits human motor cortex when applied with mostly monophasic (anterior-posterior) pulses.**

Sommer M, Ciocca M, Hannah R, Hammond P, Neef N, Paulus W, Rothwell JC (2014)

## Repetitive cTMS Specifications

Note :The multi- parameter nature of cTMS makes it difficult to describe in a table. Feel free to contact us.

|   | Pulse Type<br>(+ve pulse width,<br>-ve Pulse width) | M-Ratio | 10 Hz | 25 Hz | 50 Hz | 100 Hz | 200 Hz | 400 Hz | 800 Hz | 1000 Hz | Max Freq. (Hz)<br>at 100% output | Maximum<br>Output (%) |
|---|---|---------|-------|-------|-------|--------|--------|--------|--------|---------|----------------------------------|-----------------------|
|   |   |         |       |       |       |        |        |        |        |         |                                  |                       |
| Unidirectional train<br>Balanced Pulses | (45µs, 145µs)                                       | 0.2     | 100   | 100   | 49    | 24     | 11     | NA     | NA     | NA      | 24                               | 100                   |
|   | (60µs, 185µs)                                       | 0.22    | 100   | 82    | 40    | 19     | 9      | NA     | NA     | NA      | 20                               | 100                   |
|   | (75µs, 225µs)                                       | 0.25    | 93    | 55    | 27    | 13     | 6      | NA     | NA     | NA      | 15                               | 93                    |
|   | (11µs, 54µs)  | 0.76    | 100   | 100   | 100   | 100    | 100    | 100    | 65     | 49      | 550                              | 100                   |
|   | (20µs, 78µs)  | 0.139   | 100   | 100   | 100   | 100    | 91     | 43     | 19     | 14      | 184                              | 100                   |
|   | (30µs, 105µs)                                       | 0.1627  | 100   | 100   | 100   | 91     | 44     | 20     | NA     | NA      | 91                               | 100                   |
|   | (40µs, 131µs)                                       | 0.182   | 100   | 100   | 100   | 54     | 26     | 12     | NA     | NA      | 55                               | 100                   |
|   | (50µs, 158µs)                                       | 0.1987  | 100   | 100   | 70    | 34     | 16     | NA     | NA     | NA      | 35                               | 100                   |
|   | (60µs, 185µs)                                       | 0.2165  | 100   | 97    | 47    | 23     | 10     | NA     | NA     | NA      | 24                               | 100                   |
|   | (70µs, 212µs)                                       | 0.2365  | 97    | 66    | 32    | 16     | 7      | NA     | NA     | NA      | 17                               | 97                    |
| Bi-Directional                          | (80µs, 238µs)                                       | 0.2599  | 78    | 57    | 28    | 13     | 6      | NA     | NA     | NA      | 18                               | 85                    |
|   | (86µs, 255µs)                                       | 0.2756  | 78    | 47    | 23    | 11     | NA     | NA     | NA     | NA      | 15                               | 78                    |
|   | +ve (60µs, 100µs)                                   | 1       | 43    | 43    | 43    | 26     | 13     | 5      | NA     | NA      | 63                               | 43                    |
|   | -ve (60µs, 100µs)                                   | 1       | 43    | 43    | 43    | 21     | 10     | 4      | NA     | NA      | 50                               | 43                    |

Higher output possible with possible decay (displayed during planning step)



**Rogue Research Inc.**

6666 St-Urbain., Suite 300, Montreal (QC) Canada H2S 3H1

Telephone +1 (514) 284-3888

Toll free (North America) 866-984-3888 . Fax +1 (514) 284-6750

[www.rogue-research.com](http://www.rogue-research.com)

[info@rogue-research.com](mailto:info@rogue-research.com)

**Elevate TMS is approved for cortical and peripheral nerve stimulation for clinical use in Canada. It is not approved for clinical applications outside of Canada and is for scientific research and/or investigational use only. Operation of Elevate TMS must be done in accordance with your relevant medical device regulations.**